

microFET[®] Hand Evaluation Kit

The Hoggan Hand Evaluation Kit includes microFET digital HandGRIP and Pinch Gauge to offer handgrip and pinch strength testing in a single convenient carrying case.



MKT-008.A

microFET[®] Hand Evaluation Kit

The Hoggan Hand Evaluation Kit includes microFET digital HandGRIP and Pinch Gauge to offer handgrip and pinch strength testin in a single convenient carrying case.

microFET Pinch Features

- Load cell sensor technology providing more accurate grip and pinch strength measurements
- Easy to ready LCD display shows peak force and elapsed time
- Durable, lightweight pinch device.
- Perform tests as a standalone device, or wireless with available MicroFET or data collection software
- Selectable units of measure: pounds (lbs.), Newtons (N) or kilogram force (kgf)
- Measurement accuracy within 1%
- Auto shut-off to conserve energy
- Uses rechargeable lithium-ion battery

Hand Evaluation Kit Specs

Pinch Specs

- Measure tip, key and palmar pinch strength
- Starts recording after 0.8 lb. force to 60 lbs. in 0.1 lb. increments

HandGRIP Specs

- Can perform Maximum Grip test, 5 Position Grip test, Rapid Exchange Grip (REG).
- Low and high threshold settings:
Low – 0.8 lb. to 200 lbs. in 0.1 lb. increments
and High – 3 lbs. to 200 lbs. in 0.1 lb. increments

Your Purchase Includes

- MicroFET Pinch
- HandGRIP gauge
- Wall pack power supply (For both HandGRIP and Pinch)
- Carrying case
- User Manual
- Calibration Certificate
- 1-year standard warranty with extended warranties available
- Optional MicroFET Clinical Software and FDC FET data collection software available