

## About the microFET HandGRIP

microFET HandGRIP dynamometer measures grip forces to test and record grip strength deficits, evaluate consistency of effort, and to perform general handgrip strength assessments. microFET HandGRIP provides accurate, objective and quantifiable results for identifying areas of weakness in the wrist and hand due to injury or disease.





MKT-003.A



1.800.678.7888 | WWW.HOGGANSCIENTIFIC.COM | 1987 SOUTH 3653 WEST, SALT LAKE CITY, UTAH 84104





#### Features

- Features load cell sensor technology providing more accurate grip and pinch strength measurements
- Easy-to-read LCD display shows peak force and elapsed time
- Durable, lightweight devices
- Perform tests as a standalone device, or wireless with available MicroFET or FET data collection software
- Measurement accuracy to within 1%
- Auto shut-off to conserve energy
- Uses rechargeable lithium-ion battery

# Specifications

### **Pinch Specs**

- Measure tip, key and palmar pinch strength
- Starts recording after 0.8 lb. force to 60 lbs. in 0.1 lb. increments

#### HandGRIP Specs

(SO 13485 Guilty Madar Davies CERTFIED

(F

- Can perform Maximum Grip test, 5 Position Grip test, Rapid Exchange Grip (REG).
- Low and high threshold settings:
  Low 0.8 lb. to 200 lbs. in 0.1 lb. increments and
  High 3 lbs. to 200 lbs. in 0.1 lb. increments

## Your Purchase Includes

- MicroFET Pinch
- HandGRIP gauge
- Wall pack power supply (For both HandGRIP and Pinch)
- Carrying case
- User Manual
- Calibration Certificate
- 1-year standard warranty with extended warranties available
- Optional MicroFET Clinical Software and FDC FET data collection software available

